









100

A Single Delivery System for Right-Handed and Left-Handed Dentistry

Whether you are just starting out, considering the addition of a new operatory, or planning a multi-room clinic, Belmont Swing-Mounted Units adapt to the changing needs of your practice.

Doctor's 3 Handpiece Control with syringe & bottled water system.



A low profile Assistant's Instrument arm allows for convenient placement of evacuation valves & syringe.

X-Calibur Swing-Mounted
Delivery Systems offer you
a variety of unit configurations
including combinations of
Doctor's Control, Clesta[™]
Light, Assistant's Instruments
and Cuspidor.

These components are mounted on arms that can be re-positioned from the right side of the chair to the left side within seconds.

Versatile Equipment for the Modern Dental Operatory



Stream-lined, flat panel computer monitors give a sophisticated, high tech look to your dental practice. Belmont offers a monitor arm* that can be used alone on a chair or a version that can be coupled with a delivery system.

Select a Belmont swingmounted delivery system and match it with a Belmont X-Calibur or AcuTrac chair.





21

Develop your blueprint for the future by contacting your authorized Belmont dealer today.





Delivery Systems									
	BDS-2575	BDS-2570	BDS-2580	BDS-258D	BDS-0020	BDS-0050	BDS-0099	BDS-2530	BDS-2590
SWING MOUNTED UNITS,	Dr. Control	Dr Control &	Dr Control,	Dr Control,	Assistant	Assistant	Assistant	Clesta Light	Monitor Art
SWING MOUNTED Clesta TM LIGHT,		Clesta Light	Clesta Light, Cuspidor,	Cuspidor and Assistant	Vac Pac for Bel-20	Vac Pac for Bel-50	Vac Pac for S-99	with Chair Adapter	with Chair Adapter
SWING MOUNTED ACCESSORIES			Assistant Instruments	Instruments (No Light)	for Ber 20		AcuTrac	(only)	(only)
B Handpiece Automatic Control									
Individual Handpiece Pressure Adjustment									
Individual Handpiece Water Coolant Adjustment									
Master Handpiece Air Coolant Adjustment									
Non-retracting Water Coolant with Master On/Off									
Handpiece Flush									
4 Hole Straight Gray Asepsis Tubing with Midwest Style Connector									
1 Doctor Syringe	-								
Balance Arm with Pneumatic Brake	-								
Remote Toggles For Cup Fill & Bowl Rinse	_	_							
Swing Arm with Clesta™ light on a Curved Post			-	_					
Swing Arm for Flat Panel Monitor		_	_					_	
Factory Installed Wire Harness for Add-on Accessories									_
Factory Installed Wire Harness for Unit Mounted Touchpad			-						
Built-In City/Bottled Water System			-						
Vacuum Canister With Disposable Solids Collector Screen	_	_	_	-					
1/4' Water Outlet With Flow Control				-			_		
1/4" Air Quick Disconnect			-	_	_	_	_		
Wet/Dry Foot Control with Signal Relay			_	-					
Stainless Steel Tray with Non-Slip Pad	-			-					
	_	_	_	_					
Gravity Drain Cuspidor									
Swing Away Bowl									
Adjustable Timed Bowl Rinse									
Top Mounted Cup Fill & Bowl Rinse Buttons									
Assistant 's Instruments									
4 Position Crescent Holder									
1 Assistant Syringe, 1 HVE, 1 SE									
	1	1							
Belmont Chair Adapter (Specify)	_				D-1-20	Del 50		-	
Bel-20, Bel-50 Type					Bel-20	Bel-50			
Acutrac Type									
	_		S-99 Version of AcuTrac	S-99 Version of AcuTrac			S-99 Version of AcuTrac	_	
Equipment Utility Center	_	•						-	
Equipment Utility Center (J-Box cover & frame not shipped with units installed on	-							_	
1			of AcuTrac	of AcuTrac				Cover & Frame Only	
(J-Box cover & frame not shipped with units installed on AcuTrae Chairs w/ Long Pump Cover, Internal Umbilical)			of AcuTrac	of AcuTrac				Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options	•		of AcuTrac	of AcuTrac				Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece	•		of AcuTrac	of AcuTrac			of AcuTrac	Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side)	•		of AcuTrac	of AcuTrac				Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side) BDS0011 Additional HVE	•		of AcuTrac	of AcuTrac	•		of AcuTrac	Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side) BDS0011 Additional HVE 90-175 Additional HVE	•		of AcuTrac	of AcuTrac			of AcuTrac	Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side) BDS0011 Additional HVE 90-175 Additional HVE BDS0014 Assistant Short Arm (Replaces Telescoping Arm)			of AcuTrac	of AcuTrac			of AcuTrac	Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side) BDS0011 Additional HVE 90-175 Additional HVE BDS0014 Assistant Short Arm (Replaces Telescoping Arm) BDS0017 Delete 1 Handpiece			of AcuTrac	of AcuTrac			of AcuTrac	Cover & Frame Only	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side) BDS0011 Additional HVE 90-175 Additional HVE BDS0014 Assistant Short Arm (Replaces Telescoping Arm) BDS0017 Delete 1 Handpiece Additional Umbilical			of AcuTrac	of AcuTrac			of AcuTrac	Cover &	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side) BDS0011 Additional HVE 90-175 Additional HVE BDS0014 Assistant Short Arm (Replaces Telescoping Arm) BDS0017 Delete 1 Handpiece Additional Umbilical FPMM Flat Panel Monitor Mount			of AcuTrac	of AcuTrac			of AcuTrac	Cover & Frame Only	
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical) Available Options BDS0004 Additional Handpiece BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side)			of AcuTrac	of AcuTrac			of AcuTrac	Cover & Frame Only	



7950 Cherry Ave. Suite 103, Fontana, CA 92336 866.274.8351 • f 951.848.9700 www.ashteldental.com

Specifications, prices and models are subject to change without notice. BRH00015 REV082506